


jjinstruments.com



To review our full line
of products or
to order a catalog,
visit our website
at jjinstruments.com.

Silk™



The Premium Line
Of New Lightweight,
Superb Handling
Curettes And Scalers
With One Huge
Advantage ...

... SILK Currettes And Scalers Cost A Lot Less!

Compare our instruments with the best in the industry and you will immediately recognize we are equal in every way but one...price. Apples to apples. Our scalers and currettes to their scalers and currettes. Note the deep precision knurling for a better tactile feel. Same high-grade German 440A stainless steel tips. Same 9.5 mm diameter for easier handling. Same light weight (under 20 grams) for less hand fatigue. Same professional craftsmanship. Same high industry standards.

Of course, any manufacturer can make these claims. We guarantee and certify ours.


highest standards.

The only difference is our price. We cost a lot less.

[27 Of The Most Popular Patterns.]

In fact, our instruments are manufactured within the guidelines of the FDA quality assurance system for Good Manufacturing Practices and our factory is ISO 13485 certified. All our products bear the CE mark, emblematic of the European Community's

And if you don't agree ... if our instruments don't meet your standards, you can return them for a full refund. Decide today ... ask your dealer for *Silk* scalers and currettes.



- Beautiful, glass-bead blasted finish on center bar.
- Deep cross-cut, knurling for better non-slip control.
- Our factory is ISO 13485 certified.
- CNC machined for excellent balance, precision and consistency.
- 9.5 mm diameter for easier handling and better tactile sensitivity ... better comfort throughout the day.
- Evenly balanced for consistent control when switching from tip to tip.

- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for Good Manufacturing Practices.
- Highest grade German 440A stainless steel, razor sharp, durable tips.
- Rockwell hardened to 60 RHc.
- 27 of the most popular designs, from Graceys to McCalls to Columbias, and many European patterns.
- All instruments are guaranteed against defects in material and workmanship.
- All instruments bear the CE mark, emblematic of the European Community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... why pay more for the same high-quality instruments?

