

# ORTHODONTIC CONTENTS

---

- Bracket Placers ..... 136
- Cotton & Dressing Pliers ..... 136-137
- Band Pushers & Seaters ..... 138
- Ligature Tie Pliers ..... 139-142
- Orthodontic Instruments ..... 143
- Hollow Core Ligature Orthodontic  
Instruments ..... 143-144
- Elite Orthodontic Pliers ..... 145-150
- Orthodontic Pliers ..... 151-164

ORTHODONTIC



ORTHODONTIC



**5.25" Anterior**  
Bracket Placer

55-0008



**4.25" Posterior**  
Bracket Placer

55-0009



**College Plier**  
Serrated Tips

08-200



**College Plier**  
Serrated, Narrow Tips

08-205



**College Plier**  
Self-Locking, Serrated

08-210



**College Plier**  
Self-Locking, Serrated,  
Narrow Tips

08-215



**College Plier**  
Self-Locking, Grooved

08-220



**Markley Plier**  
6" Self-Locking, 1 Grooved Tip,  
1 Serrated Tip

08-250



**Meriam Plier**  
6.25" Serrated

08-255



**Perry Plier**  
5" Serrated

08-262

ORTHODONTIC

# Band Pushers & Seaters

ORTHODONTIC



Band Pusher/  
Ligature Director

07-815



**#300S**  
Band Pusher/Scaler

07-820



**#300**  
Band Pusher

07-825



**#1400**  
Band Seater  
W/ Tin Insert

07-830



**6" Autoclavable (290°)**  
Plastic Band Seater  
W/ Stainless Steel Insert

07-835



Actual Size

**5.5" Regular Tip**  
Mathieu Stainless Steel

25-7650



**5.5" Regular Tip**  
Mathieu Gator Grip™  
W/ Carbide Inserts

65-7650



**5.5" Delicate Tip**  
Mathieu Gator Grip™  
W/ Carbide Inserts

65-7655



**5.5" Delicate Tip**  
Mathieu  
Stainless Steel

25-7655



**5.5" Elastomeric**  
Mathieu

25-7660



Actual Size

**5.5" W/ Notched Tip**  
Mathieu

25-7665



**5.5" W/ Hook**  
Mathieu

25-7670

ORTHODONTIC

# Ligature Tie Pliers



**4.75" Delicate Tip**  
Mathieu Stainless Steel  
25-7680



**4.75" Delicate Tip**  
Mathieu Gator Grip™  
W/ Carbide Inserts  
65-7680



**4.75" Elastomeric**  
Mathieu  
25-7675



**4.75" W/ Hook**  
Mathieu  
25-7685



**5" STR**  
Mosquito FCP  
24-1050



**5" CVD**  
Mosquito FCP  
24-1055




**5" STR W/ Hook**  
Mosquito FCP  
24-1052



**5" CVD W/ Hook**  
Mosquito FCP  
24-1057

ORTHODONTIC



6"  
Crile-Wood N/H  
25-6420



6" Carbide Ester Grip™  
Crile-Wood N/H  
65-6420

ORTHODONTIC

# Ligature Tie Pliers



Actual Size



**5" Elastomeric**

Boynton N/H

25-7550



**5.5" With Notch**

Mathieu N/H

25-7665



**5" Carbide**

Boynton N/H

65-7550

ORTHODONTIC



J&J Instruments is proud to offer our Hollow Core™ orthodontic instruments. The handles measure a full 3/8" in diameter and the hollow design creates an instrument that is user friendly. They fit in the hand easily and securely and so promote tactile sensitivity while reducing hand fatigue.



**Ligature**  
Director SE  
07-810



**Ligature**  
Director DE  
07-812



**Ligature**  
Director/Scaler DE  
07-817



**Hollow Core™**  
Ligature Director/Scaler  
77-817




**Hollow Core™**  
Ligature Director/Director  
77-819

ORTHODONTIC

# Hollow Core™ Ligature Directors & Band Pushers

ORTHODONTIC



**Hollow Core™**  
Ligature Director SE

77-810

**#H5 Hollow Core™**  
Band Pusher/Scaler

77-822

**DE Hollow Core™**  
Ligature Director DE

77-812

**Hollow Core™**  
Band Push/Lig Director

77-815

**#300S Hollow Core™**  
Band Pusher/Scaler

77-820

The Elite Plier line is manufactured from a high strength European steel. Each plier comes with its own Allen wrench, and the special screw joints can be adjusted in the office so the "feel" will always be the way you want it. The fluted handle design promotes comfort while assuring a positive grip.



**#16 Reg. Head**  
Distal End Cutter

29-016




**#16S Small Head**  
Distal End Cutter

29-017



Distal End Cutter with safety hold and inserted tool steel tips. Cuts hard, round and rectangular arch wires. Holds loose distal end wire.




**#1002**  
Pin/Lig Cutter

29-1002



Pin & Ligature Cutter with inserted tool steel tips for close cutting of pins and ligatures.

ORTHODONTIC



**#101**  
Force Module Separator

29-101

Force Module Separator to easily hold, spread, and place elastics. Adjustable screw for control of width opening.



**#110**  
How Plier

29-110

How Plier. Straight utility plier with serrated pads and slip-safe beaks.



**#139 Short**  
Bird Beak Plier

29-139



**#141 Long**  
Bird Beak Plier

29-141



**#158**  
Weingart Plier

29-158



**#200**  
Clasp Adjusting Plier

29-200



Bird Beak. Snub nose. One round cone beak, one pyramid beak. For working round wire.



Weingart Utility. Utility plier with offset tips and serrated pads.



Clasp Adjusting. Three jaw plier. Straight aligned tips. For delicate, consistent bending. Used in the closing and adjusting of clasps.



**#347 Short**  
Posterior Band Remover

29-347



**#347 Long**  
Posterior Band Remover

29-349


Posterior Band Remover. Plier nose on tip without pad has been shortened. Replacement tips available.

Posterior Band Remover. Standard posterior band remover. The plier nose on the non-padded tip is longer than the 29-347. Replacement tips available.



**#348**  
Direct Bracket Remover  
29-348

Direct Bracket Remover. The plier nose is thinner than the nose on the 29-347. Removes wedge between both edges of tooth surface. Replacement tips available.



**#352**  
Omega Loop Forming Plier  
29-352  
**4/PK**  
Band Remover Replacement Tips  
29-998

Omega Loop Forming Plier. Fluted handle for sure grip, adjustable screw lock, with included wrench.

ORTHODONTIC



**#7001 Elite**  
Wire Cutter

29-7001




**Elite**  
Screw & Wrench Kit

29-999



Our orthodontic pliers are perfect for the General Dentist and Cosmetic Dentist who require instruments to make minor adjustments to removable prosthetics such as partials and retainers. The wide variety of patterns we offer all feature a smooth matte finish, uniformly beveled edges, and precisely formed beaks. The J&J Instrument assurance of excellent quality and exceptional value makes our orthodontic pliers a sensible addition to every office's instrument inventory.



**#150**

Pin & Ligature Cutter

09-150



Pin and Ligature Cutter. For close cutting of pins and ligatures.

**#150S Small**

Pin & Ligature Cutter

09-151



Pin and Ligature Cutter Small Handle. For close cutting of pins and ligatures. For small hands.




**#150 Carbide**

Pin & Ligature Cutter

69-150



Pin and Ligature Cutter Carbide Tips. For close cutting of pins and ligatures with tungsten-carbide inserted tips.



**#151S Carbide**  
Pin & Ligature Cutter

69-151



**#16**  
Distal End Cutter

09-016



Distal End Cutter. With safety hold and inserted steel tips. Cuts hard round and rectangular arch wire. Holds loose distal end wire.



Diagonal Hard Wire Cutter

09-177



**#16S Small**  
Distal End Cutter

09-017

Small head with safety hold and inserted tool steel tips. Cuts hard round and rectangular arch wire. Holds loose distal end wire. Smaller head than 09-016.



Arrow Head Forming Plier

09-000

V-form wire contouring plier with double notched tip.



**#12**  
Young Wire-Bending Plier

09-012

For orthodontic or lab use. Makes various bends and loops for wire bending.

# Orthodontic Pliers

ORTHODONTIC



**#53**  
Wire-Bending Plier

09-053

One concave beak and one round beak. Used to bend clasps and wires.



**#70**  
Side Cutting Plier

09-070

Cuts round hardwire. Cutting angle 15° to improve vision and accessibility.



**#100**  
Lingual Bar & Wire Bending Plier

09-100

To bend and shape lingual bars and wires.



**#101**  
Wire Bending 3 Groove Plier

09-101

For bending wire.



**#102**  
Cone & Socket Plier

09-102

To tighten and remove cone socket points.



**#103**  
Force Module Separator

09-103

To easily hold, spread, and place elastics. It is lightweight. No adjusting screw.



**#104**  
Flat Nose Plier

09-104

For clasp and wire bending.



**#107**  
Round Nose Plier

09-107

Smooth, tapered beaks. For shaping crown and bands. For wire bending.



**#110**  
How Straight Plier

09-110

Straight utility plier, serrated pads, slip safe beaks.



**#111**  
How Curved Plier  
09-111

45° utility plier, serrated pads.



**#112**  
Abell Contouring Plier  
09-112

Ball and socket beaks. Used for crown and band contouring.



**#114**  
Johnson Contouring Plier  
09-114

For contouring bands and crowns.

ORTHODONTIC



**#115**  
Reynolds Contouring Plier

09-115

Falcon beaks. For contouring bands and crowns.



**#118**  
Peeso Collar Plier

09-118

One flat beak and one convex beak. Used to shape and fit collars or cap crowns.



**#120**  
Robinson Collar Plier

09-120

For band shaping and bending.





Serrated

Smooth

**#121 Serrated**

Lab/Office Plier

09-121

For wire bending and general use. Serrated beaks.

**#122 Smooth**

Lab/Office Plier

09-122

For wire bending and general use. Smooth beaks.



**#134**

McKellop Clasp Plier

09-134

For clasp bending.



**#136**

Goslee Contour Plier

09-136

One beak smooth, one beak serrated. General use plier.



**#137**  
Gordon Contouring Plier

09-137

For general contouring.



**#139 Short**  
Bird Beak Plier

09-139

Snub nose. One round cone beak, one pyramid beak. For working round wire.



**#140 Long**  
Bird Beak Plier

09-140

Snub nose. One round cone beak, one pyramid beak. For working round wire. Longer taper beaks than the 09-139. For lighter wire.



**#142A**  
Ribbon Arch Forming Plier

09-142

To hold and bend ribbon arch wire.



**#156**  
Fischer Flat Nose Plier

09-156

Flat nose lip. For general use and for ligating. Serrated beaks.



**#160**  
Anterior Band Removing Plier

09-160

For removing anterior bands and brackets.

ORTHODONTIC

## Orthodontic Pliers



**#200**  
Clasp Adjusting Plier

09-200

Three prong plier. Straight aligned tips. For delicate, consistent bending. Used in the closing and adjusting of clasps.



**#347**  
Oliver Band Remover

09-347

To remove anterior and posterior bands.



**#301**  
Triplex Contouring Plier

09-301


Three in one: Johnson, Reynolds, and band stretching.



Tweed Loop Forming Plier

09-350

For making bull and helical closing loops in square or rectangular wires.



#355

Jarabak Plier

09-355

For bending lighter wire. For a firm grip. Flat tip has serrations.




#417

Crown & Shell Plier

09-417

Crown and shell crimping.



**#858**  
Dela-Rosa Plier

09-858

Forms inner bows and loops.



**#861**  
Weingart Utility Plier

09-861

Utility plier. Offset tips. Serrated, pointed pads hold wire firmly at a convenient working angle.