


RESTORATIVE CONTENTS

<i>Silk</i> Composite Instruments	110
Plastic Filling Instruments	111-112
Excavators	113-114
Cavity Liners	115
Amalgam Carriers	116-117
Carvers	118-120
Cleoid Discoids	120
Pluggers / Condensers / Endodontic Spreaders	121-123
Burnishers	124-126
Articulating Paper Forceps	127
Cement Spatulas	128-129
Wax Spatulas	130
Wax Instruments	131
Calipers and Gauges	132
Crown Removers	133
Matrix Retainers	134
Rubber Dam Instruments	134-135



The *Silk* Composite Instruments follow the designs of the most popular restorative instrument patterns with some notable improvements. The tips are a bit thinner and more delicate. This enhances the user's ability to sculpt and facilitates interproximal manipulation. The tips are given an ultra high polish finish which inhibits sticking and pullback of composite materials. These improved tips are then fitted into our new lightweight *Silk* handle to promote comfort and precision in every procedure.

Your restorations will be as smooth as, well, *Silk*.

RESTORATIVE



#3
P.F.I.
95-030



#2
P.F.I.
95-020



#21B
Acorn Burnisher
95-600



#3
Woodson P.F.I.
95-103



#3S
Hollenback Carver
95-230



#3
P.F.I.
08-030



#4/5
Gregg P.F.I.
08-000



#2
Ladmore P.F.I.
08-002



#3
Ladmore P.F.I.
08-003



#1
Woodson P.F.I.
08-101



#2
Woodson P.F.I.
08-102



#3
Woodson P.F.I.
08-103

RESTORATIVE

Plastic Filling Instruments



#13
P.F.I.
08-060



#1
P.F.I.
08-010



#2
P.F.I.
08-020



#4
P.F.I.
08-040



#5
P.F.I.
08-050



#11
P.F.I.
08-055

RESTORATIVE



#5
Excavator
04-705



#6
Excavator
04-706



#14
Excavator
04-714



#17
Excavator
04-717



#18
Excavator
04-718




#19
Excavator
04-719



#20
Excavator
04-720



#22
Excavator
04-722



#24
Excavator
04-724


RESTORATIVE

Excavators



#38/39
Excavator

04-738



#61/62
Excavator

04-761



#31LS Spoon
Excavator

04-730



#31L
Excavator

04-731



#32L
Excavator

04-732



#33L
Excavator

04-733



#2 Glick
Excavator

04-765

RESTORATIVE



4.25" Regular Tip
Cavity Liner

02-540



5.5" Long Tip
Cavity Liner

02-545



6" DE Reg/Long
Cavity Liner

02-550

RESTORATIVE

Amalgam Carriers

RESTORATIVE



DE Reg/Lg
Amalgam Carrier

01-030



DE Reg/Reg
Amalgam Carrier

01-040



DE Reg/Jum
Amalgam Carrier

01-050



DE Lg/Lg
Amalgam Carrier

01-060



DE Jum/Jum
Amalgam Carrier

01-070



Regular Large Jumbo Gun Type

Amalgam Carrier tips shown actual size.



SE Regular
Amalgam Carrier

01-300



SE Large
Amalgam Carrier

01-412



SE Jumbo
Amalgam Carrier

01-423



SE Gun Type
Amalgam Carrier

01-500

RESTORATIVE



Regular Large Jumbo Gun Type

Amalgam Carrier tips shown actual size.

RESTORATIVE



Knife Tip Only
Akers Carver

03-193

Handle Only
Akers Carver

03-191

Complete
Akers Carver

03-190

Spoon Tip Only
Akers Carver

03-192



Le Cron Carver

03-300



Reeves Carver

03-320



Roach Carver

03-340



Vehe Carver

03-360




#3S (1/2)
Hollenback Carver
03-230



Frahm Carver
03-195



#3
Hollenback Carver
03-200



Martin Carver
03-310



#7
Levy Carver
03-330



Interproximal Carver
03-335



Shoshan-A Carver
03-350



#1
Walls Carver
03-410

RESTORATIVE

Carvers

Cleoid Discoids



#3
Walls Carver
03-430



#1
Ward Carver
03-450



#1S (1/2)
Ward Carver
03-460



#2
Ward Carver
03-470



#89/92
Cleoid Discoid
04-689



#3/6 Large
Cleoid Discoid
04-603



#WACD
Cleoid Discoid
04-610



#90/93
Cleoid Discoid
04-690

RESTORATIVE



#1/2 Serrated
Behrman Amalgam
Plugger

01-220



#0/1 Serrated
Black Plugger

01-210



#0/1 Smooth
Black Plugger

01-215



#1/2 Smooth
Black Plugger

01-225



#2/3 Serrated
Black Plugger

01-240



#2/3 Smooth
Black Plugger

01-245



#4/5 Serrated
Black Plugger

01-246




#4/5 Smooth
Black Plugger

01-247



#9/10
Back Action Plugger

01-200



Black Tip Sizes

RESTORATIVE



#11/12
Back Action Plugger

01-205



Smith Special
Amalgam Plugger

08-090




#1
Behrman Amalgam Plugger

01-010



#1 Smooth
Hollenback Condenser

01-250



#1 Serrated
Hollenback Condenser

01-255




#2 Smooth
Hollenback Condenser

01-260




#2 Serrated
Hollenback Condenser

01-265



#1T JB
Tanner DE Condenser

01-410



#2T
Tanner DE Condenser

01-420



#3T
Tanner DE Condenser

01-430



#5T
Tanner Carver DE

01-450



#D11
Root Canal
Spreader

04-770

#D11T
Root Canal
Spreader

04-772

#MA57
Root Canal
Spreader

04-774

#3
Root Canal
Spreader

04-776

RESTORATIVE

Burnishers

RESTORATIVE



#28/29 Smooth
Burnisher DE
02-500



#2/3
Burnisher DE
02-501



#2/26
Burnisher DE
02-502



#2/34
Burnisher DE
02-503



#25/26 (13A)
Burnisher DE
02-504




#26/27S
Burnisher DE
02-505



#29/34
Burnisher DE
02-506



#26/29
Burnisher DE
02-508



#27/29
Burnisher DE
02-510



#21B
Acorn Burnisher DE
02-600



#21
Acorn Burnisher DE
02-522



#21C
Acorn Burnisher DE
02-601



#21D
Acorn Burnisher DE
02-602



#1
Ball Burnisher DE
02-610



#1
U/Michigan Burnisher
02-710



#2
U/Michigan Burnisher
02-720

RESTORATIVE

Burnishers

RESTORATIVE



#21
Acorn Burnisher SE

02-521



#28
Burnisher SE

02-528



#29
Burnisher SE

02-529



#33
Burnisher SE

02-533




#34
Burnisher SE

02-534


Articulating Paper Forceps

Precision crafted for years of dependable service. Manufactured from a high strength type 420A superior quality European stainless steel, our articulating paper forceps will withstand repeated autoclaving. Longitudinal grooves run the length of the high tension jaws. Cross hatch serrations at the tip assure a tight no-slip grip on the articulating paper. The tongue and grooved feature of the 08-305 is designed to hold even the thinnest Mylar® film in place.




Articulating Paper Forceps Straight Serrated

08-300



Articulating Paper Forceps Straight Grooved

08-305



Articulating Paper Forceps Curved Serrated

08-310

Cement Spatulas

RESTORATIVE



#5
Cement Spatula DE

13-105



#7
Cement Spatula DE

13-107



#26
Cement Spatula DE

13-126



#11S (313)
Cement Spatula

13-110



#22
Cement Spatula

13-122



#24
Cement Spatula

13-124



#24A
Cement Spatula

13-125



#324 Rigid
Cement Spatula

13-324



#336
Cement Spatula

13-336

RESTORATIVE

Wax Spatulas

RESTORATIVE



#5
Wax Spatula

13-250



#7
Wax Spatula

13-270



#7A St. Louis
Wax Spatula

13-275



#7.5 Baltimore
Wax Spatula

13-277



#31
Wax Spatula

13-310



#45
Wax Spatula

13-450



#1
PK Thomas Wax
13-301



#2
PK Thomas Wax
13-302



#3
PK Thomas Wax
13-303



#4
PK Thomas Wax
13-304



#5
PK Thomas Wax
13-305

Set has 1 each of all five
PK Thomas Instruments #1-5.

#1-5 (Set of 5)
PK Thomas Wax
13-306




#2S
Gold Knife U.S.C.
08-080



#2
Gold Knife U.S.C.
08-070

RESTORATIVE

Calipers and Gauges



Boley Millimeter Gauge with Lock

06-800



Iwansson Spring Caliper Blunt Tip for Wax

06-840



Wedge Shaped mm Gauge

06-850



Iwansson Spring Caliper Sharp Tip for Metal

06-841

RESTORATIVE



W/ A, B, and C Tips
Morrell Crown Remover

07-898



#A Only
Morrell Tip

07-895



#B Only
Morrell Tip

07-896



#C Only
Morrell Tip

07-897



#255
Baade Band Removing Plier

09-255

RESTORATIVE


Matrix Retainers

Rubber Dam Instruments



Tofflemire Universal
Child

07-295



Tofflemire Universal
Adult

07-300



Ainsworth Rubber
Dam Punch *

11-500



Tofflemire Contra
Angle Child

07-305



Tofflemire Contra
Angle Adult

07-310



IV-Type Rubber
Dam Punch *

11-550

* Not actual size

RESTORATIVE



Brewer Rubber Dam
Clamp FCP *

11-525




U/Wash Rubber Dam
Clamp FCP *

11-595



IV-Type Rubber
Dam Clamp FCP *

11-575



Rubber Dam
Frame 5" *

11-585

RESTORATIVE

* Not actual size